college park
TECHNOLOGY for the HUMAN RACE

technical instructions

التعليمات الفنية · Teknisk vejledning · Technische Anleitung
Τεχνικές οδηγίες · Instrucciones técnicas · Tekniset ohjeet
Instructions techniques · הוראות טכניות · Istruzioni tecniche
Technische instructies · Tekniske instruksjoner · Instrukcje techniczne
Instruções técnicas · Instruções técnicas · Технические инструкции ·
Technické pokyny · Tekniska anvisningar · Teknik Talimatlar · 技术说明

FIGURE 1




FIGURE 2


FIGURE 3



FIGURE 4

13-15 cm: 0.23 in (5.8 mm)
16-18 cm: 0.25 in (6.4 mm)
19-22 cm: 0.38 in (9.5 mm)

Package Contents

- (1) Scout Foot
(1) CPI Sock
(1) Sealing Boot (16-22 cm upon request)
- (1) Foot Shell
(1) Scout Wedge Kit

Tools Required

- (1) 3 mm Hex Key

Tools Recommended

- (1) Foot Horn

This diagram (Figure 1) is to help familiarize you with the unique parts of the Scout. These parts are referenced in the instructions and used when speaking with a technical service representative.

Key Components (Figure 1)

Foot

- A. Housing
E. Rounded Heel Wedge (optional)
- B. 3mm Clamp Screw - Torque 5 N·m (44 in-lbs)
F. Heel Spring
- C. Toe Spring
• CPI Sock (not shown)
- D. Permanent Wedge
• Foot Shell (not shown)

Threaded Adapter Options (additional purchases)

- G. Pyramid Adapter
K. Lamination Adapter - Medium
- H. Pyramid Receiver
L. Lamination Adapter - Large
- I. Exo Receiver
- J. Lamination Adapter - Small

Exoskeletal Mounting (additional purchases)

- M. Exo Pyramid Tool Kit (reusable)
- N. Exo Block Kit
- O. Exo Mounting Bolt (included w/exo block kit) - Torque 47 N·m (35 ft-lbs)

PRODUCT DESCRIPTION

This prosthetic foot device is constructed with a housing, two composite springs, and heel wedge. The toe spring is secured to the housing and heel spring with fasteners.

INTENDED USE

The Scout is a prosthetic foot designed to replace one or more functions of the biologic human foot.

INDICATIONS:

- Lower limb amputations

CONTRAINDICATIONS:

- None known

TECHNICAL SPECIFICATIONS

FOOT SIZE	WEIGHT LIMIT		
13-15 cm	55 lbs / 25 kg		
16-18 cm	110 lbs / 45 kg		
19-22 cm	132 lbs / 60 kg		
BUILD HEIGHT	13-15 CM	16-18 CM	19-22 CM
Pyramid	1.80 in / 4.6 cm	1.94 in / 4.9 cm	2.08 in / 5.3 cm
Receiver	2.0 in / 5.3 cm	2.24 in / 5.7 cm	2.38 in / 6.0 cm
Exo Block	Not Available	4.09 in / 10.4 cm	4.25 in / 10.8 cm
Lam Adapter	1.85 in / 4.7 cm	1.99 in / 5.1 cm	2.13 in / 5.4 cm

GAIT MATCHING® GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

FIRMNESS CATEGORY CHART

WEIGHT LBS	0-55	56-74	75-100	101-132
WEIGHT KG	0-25	26-33	34-45	46-60
SIZE CM	13-15			
LOW-HIGH IMPACT	1	N/A	N/A	N/A
SIZE CM	16-18			
LOW-HIGH IMPACT	1	2	3	N/A
SIZE CM	19-22			
LOW-HIGH IMPACT	1	2	3	4

ASSEMBLY AND DISASSEMBLY (FOR SOCK REPLACEMENT)

Use the Foot Horn to don and doff the foot shell. Remove the CPI sock and replace as needed. Any further disassembly or modification of components will void the warranty.

INSTALLING THREADED ADAPTER

1. Screw threaded adapter or threaded receiver completely into housing (Figure 2).
2. To set internal-external rotation, turn counterclockwise no more than 90 degrees (1/4 turn) from starting position to desired toe out angle (Figure 3). Caution: To prevent foot separation from prosthesis, do not exceed 90 degrees (1/4 turn).
3. Torque 3mm clamp screw 5 N·m (44 in·lbs).

ENDOSKELETAL MOUNTING

Use only high quality proximal endoskeletal components.

EXOSKELETAL MOUNTING

Refer to Scout Exo Block Kit Instruction Sheet for detailed instructions.

1. Attach threaded exo receiver to exo block with anti-rotation pin(s) oriented in desired position. Apply Loctite® 242 to mounting bolt. Torque to 47 N·m (35 ft·lbs). To skip alignment and lamination, go to Step 9.
2. If using the alignable exo option, attach exo pyramid tool to exo block with four 5 mm screws. Torque to 15 N·m (11 ft·lbs).
3. Attach 22 mm endo components to the exo pyramid tool and temporarily mount the socket.
4. Screw threaded exo receiver completely into housing. Torque clamp screw to 5 N·m (44 in·lbs). Don CPI sock and foot shell, then perform dynamic alignment.
5. Once alignment is complete, remove foot from exo block.
6. Mount aligned prosthesis in transfer jig. Lock socket and exo block in position.
7. Remove endo components and exo pyramid tool.
8. Use desired method to span exo block to the socket and remove from jig. Shape and laminate to desired finish. Do not remove foam from the top of the exo block. Note: Roughening the exo block foam by using 80-120 grit sandpaper will greatly improve adhesion to both the lamination and foam.
9. Re-assemble threaded exo receiver to exo block. Apply Loctite® 242 to the exo mounting bolt and torque to 47 N·m (35 ft·lbs). Re-assemble foot, don CPI sock and foot shell.

For exo growth plate assembly, refer to the Truper Growth Plate Kit Instruction Sheet (included with growth plates).

LAMINATION ADAPTER MOUNTING

Refer to Scout Threaded Lamination Adapter Fabrication Instructions for detailed instructions.

1. Temporarily attach the foot to the threaded adapter and mark the position (anterior or posterior). This is to ensure the adapter is not in the wrong orientation once fabrication is complete. Remove the foot when finished.
2. Score the surface of the lamination adapter for better adhesion.

SINGLE STAGE LAMINATION:

- Apply inner PVA bag over model.
- Add double layer of perlon stockinet, with distal carbon or nyglass reinforcement.
- Tie unidirectional carbon strips or tape through the prong holes. Secure in AP/ML position.

TWO STAGE LAMINATION:

- All surfaces must be roughened prior to applying the second lamination.
 - Transfer alignment and fix the lamination adapter to the inner socket layer with bonding agent (i.e. resin mixture, putty, adhesive, etc).
3. Apply o-ring to base of lamination adapter threads to prevent contact with resin.
 4. Apply mold release to threads of lamination cap and tighten onto adapter.
 5. If necessary, apply nyglass or carbon layers under adapter ears for reinforcement.
 6. Add double layer of perlon stockinet.
 7. Proceed with layup materials appropriate for the individual.
 8. Tie layup material into distal groove.
 9. Apply outer PVA bag and laminate socket. Ensure that resin saturates all layup material.
 10. Once cured, remove excess resin and PVA bag. Then remove lamination cap and o-ring. Note: Ensure that no resin remains on the threads.
 11. Screw threaded adapter completely into housing. Torque clamp screw to 5 N·m (44 in·lbs). Don CPI sock and foot shell.

For lamination adapter growth plate assembly, refer to the Scout Threaded Lamination Adapter Growth Plate Kit Instruction Sheet (included with growth plates).

STATIC ALIGNMENT

For optimal function, balance the patient's weight evenly between the heel and toe (Figure 4).

- Each size group has a different heel rise: **13-15 cm:** 0.23 in (5.8 mm), **16-18 cm:** 0.25 in (6.4 mm), **19-22 cm:** 0.38 in (9.5 mm)
- The load line divides the foot at 1/3 heel lever and 2/3 toe lever.

DYNAMIC ADJUSTMENTS

DESIRED RESULT	ALIGNMENT CHANGE	COMPONENT CHANGE
Firmer Toe Response	Plantarflex the Scout or move load line posterior	---
Softer Toe Response	Dorsiflex the Scout or move load line anterior	---
Firmer Heel Response	Dorsiflex the Scout or move load line anterior	Add heel wedge
Softer Heel Response	Plantarflex the Scout or move load line posterior	Remove heel wedge

**To install or remove heel wedge, refer to the Scout Heel Wedge Kit Instructions.*

USE IN WATER

The Scout has been approved for use in fresh water.

- Use caution when walking on wet surfaces.
- After the foot encounters moisture, wipe it dry using a lint-free cloth.

WARNINGS

WARNING

- Do not expose this product to corrosive materials, salt water or pH extremes.
- Any further disassembly or modification of components will void the warranty.

Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

RESIDUAL RISK STATEMENT

NOTICE OF RESIDUAL RISK

- During fitting process, ensure that CPI sock does not become pinched between foot and endoskeletal componentry.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight and/or impact level may require more frequent inspections. Soft component wear depends on the patient weight, impact level and environment. We recommend you inspect the following applicable parts for excessive wear and fatigue at each warranty inspection and replace as needed.

- Composites and Adapters
- Foot Shell
- CPI Sock
- Heel Wedge (if applied)

WARRANTY INSPECTION SCHEDULE FOR THE SCOUT: SIX MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc. and the competent authority of the Member State in which the prosthetist and/or patient is established.

**"Serious incident" is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient's, user's, or other person's state of health, (c) a serious public health threat.*

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles. Depending on patient activity this may correspond to 2-3 years of use.

ISO 10328 - LABEL

FOOT SIZE	WEIGHT LIMIT (KG)	LABEL TEXT
13-15 cm	25	n/a*
16-18 cm	45	n/a*
19-22 cm	60	ISO 10328-P3-60 kg

*Actual test loads reflect maximum body mass

ISO 10328 - "P" - "m"kg*

**) Body mass limit not to be exceeded!
For specific conditions and limitations of use, see intended use section of manufacturer's written instructions.*

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©2023 College Park Industries, Inc. All rights reserved.