

Beställningsedel College Park *Accent*

CUSTOMER #	P.O.#	DATE	
PROSTHETIST INFORMATION			
BILLING		SHIPPING (LEAVE BLANK IF SAME AS BILLING)	
FACILITY/ATTN:		FACILITY/ATTN:	
ADDRESS		ADDRESS	
CITY	STATE/PROV	CITY	STATE/PROV
COUNTRY	POST	COUNTRY	POST
PHONE	FAX	PHONE	FAX
CARRIER* <input type="checkbox"/> UPS <input type="checkbox"/> OTHER		DATE REQUIRED	TIME
		EMAIL	
PATIENT INFORMATION			
PROSTHETIST NAME		REQUISITIONER	
PATIENT ID			
NOTES			

Art.nr. Right: (IP) 761600-1XX / (DP) 761700-1XX

Transtibial Transfemoral

Art.nr. Left: (IP) 761600-0XX / (DP) 761700-0XX

Transtibial Transfemoral



ADJUSTABLE HEEL

PART ID	SIDE	SHELL COLOR	MOUNTING	BUTTON PLACEMENT	SIZE	WIDTH	FIRMNESS CATEGORY
AC	L						G
AC	R						G

Additional Accessories:
 Shelltread

Caucasian	C
Brown	B

Integrated Pyramid	IP
Dynamic Pylon + Receiver	DF
Dynamic Pylon + Adapter	DM
Dynamic Pylon + 4-hole	D4

Medial	M
Lateral	L

21-28cm	Narrow	N
	Wide	W

College Park suggests the Wide option for male patients.

1-3

WEIGHT KG	0-54	55-63	64-72	73-81	82-100	101-113
SIZE CM	21-28					25-28
LOW IMPACT	1	1	2	2	3	3
MODERATE IMPACT	1	2	2	3	3	3

Erimed International
Stensättravägen 9A
127 39 Skärholmen, Sverige
Tel: 08- 449 56 50,
order@erimed.se www.erimed.se