

Beställningssedel College Park *Soleus*

CUSTOMER #	P.O.#	DATE
PROSTHETIST INFORMATION		
BILLING		SHIPPING (LEAVE BLANK IF SAME AS BILLING)
FACILITY/ATTN:		FACILITY/ATTN:
ADDRESS		ADDRESS
CITY	STATE/PROV	CITY
COUNTRY	POST	COUNTRY
PHONE	FAX	PHONE
CARRIER* ☐ UPS ☐ OTHER		DATE REQUIRED
		TIME
		EMAIL
PATIENT INFORMATION		
PROSTHETIST NAME		REQUISITIONER
PATIENT ID		
NOTES		

Art.nr. Right: 761900-1XX

☐ Transtibial ☐ Transfemoral

Art.nr. Left: 761900-0XX

☐ Transtibial ☐ Transfemoral




M3

PART ID	SIDE	SHELL COLOR	HOUSING COLOR	SIZE	WIDTH	FIRMNESS CATEGORY
SO	L					G
SO	R					G

Caucasian	C
Tan	T
Brown	B
Jet Black	J

Perfect Pink	P
Midnight Black	M
Silver	S
Gold	G

21-30cm	
Narrow	N
Wide	W
24-30cm only	

11-76

College Park suggests the Wide option for male patients.

See next page

Additional Accessories:

- ☐ Shelltread
- ☐ Rotatable Adapter
- ☐ Rotatable Receiver
- ☐ Standard Adapter - Titanium
- ☐ Standard Receiver - Titanium
- ☐ Standard Adapter - Stainless Steel
- ☐ Standard Receiver - Stainless Steel

SOLEUS FIRMNESS CATEGORIES

WEIGHT KG	0-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125
SIZE CM	21-24								
MODERATE IMPACT	11	12	33	34	54	55	55	N/A	N/A
HIGH IMPACT	12	33	34	54	55	55	N/A	N/A	N/A
SIZE CM	25-26								
MODERATE IMPACT	12	12	32	33	54	55	75	76	N/A
HIGH IMPACT	12	32	33	54	55	75	76	76	N/A
SIZE CM	27-30								
MODERATE IMPACT	11	11	12	32	33	34	54	55	55
HIGH IMPACT	12	12	32	33	34	54	55	55	55

CUSTOM ENGRAVINGS Check one box under design. Allow two extra days for processing. No additional charge. Soleus only.

NO ENGRAVING	I AM FEARLESS	I AM STRONG	I AM DYNAMIC	HEART AND SOLEUS	WHOOOPS A DAISIES	SUPER STARS	LIGHTNING ON YOUR FEET	HOME OF THE BRAVE	FLAME ON
									