

Beställningsedel College Park *Trustep*

CUSTOMER #	P.O.#	DATE
PROSTHETIST INFORMATION		
BILLING		SHIPPING (LEAVE BLANK IF SAME AS BILLING)
FACILITY/ATTN:		FACILITY/ATTN:
ADDRESS		ADDRESS
CITY	STATE/PROV	CITY
COUNTRY	POST	COUNTRY
PHONE	FAX	PHONE
CARRIER* <input type="checkbox"/> UPS <input type="checkbox"/> OTHER		DATE REQUIRED
		TIME
		EMAIL
PATIENT INFORMATION		
PROSTHETIST NAME		REQUISITIONER
PATIENT ID		
NOTES		

Art.nr. Right: 761100-1XX

Transtibial Transfemoral

Art.nr. Left: 761100-0XX

Transtibial Transfemoral



trustep **M3**

PART ID	SIDE	SHELL COLOR	MOUNTING	SIZE	WIDTH	FIRMNESS CATEGORY
TS	L					G
TS	R					G

Additional Accessories:
 Shelltread
 Exo Block
 Exo Alignment Tool
 Endo Sealing Boot

Caucasian	C	Endo	EN	22-31cm	Narrow	N	1-9
Tan	T	Exo	ALX		Wide	W	
Brown	B				24-31cm only		
Jet Black	J				College Park suggests the Wide option for male patients.		

WEIGHT KG	0-33	34-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125	126-136	137-150	151-160
SIZE CM	22-31											EXO ONLY	
LOW IMPACT	1	2	3	4	4	5	6	6	7	7	8	8	9
MODERATE IMPACT	1	2	3	4	5	6	6	7	7	7	8	9	9
HIGH IMPACT	2	3	4	5	6	6	7	7	8	8	8	9	9

Erimed International
 Stensättravägen 9A
 127 39 Skärholmen
 Tel: 08- 449 56 50,
order@erimed.se www.erimed.se