

Beställningsedel College Park *Velocity*


CUSTOMER #	P.O.#	DATE	
PROSTHETIST INFORMATION			
BILLING		SHIPPING (LEAVE BLANK IF SAME AS BILLING)	
FACILITY/ATTN:		FACILITY/ATTN:	
ADDRESS		ADDRESS	
CITY	STATE/PROV	CITY	STATE/PROV
COUNTRY	POST	COUNTRY	POST
PHONE	FAX	PHONE	FAX
CARRIER* <input type="checkbox"/> UPS <input type="checkbox"/> OTHER		DATE REQUIRED	TIME
		EMAIL	
PATIENT INFORMATION			
PROSTHETIST NAME		REQUISITIONER	
PATIENT ID			
NOTES			

Art.nr. Right: 761910-1XX

Transtibial Transfemoral

Art.nr. Left: 761910-0XX

Transtibial Transfemoral



velocity **M3**

PART ID	SIDE	SHELL COLOR	MOUNTING	SIZE	WIDTH	FIRMNESS CATEGORY
VL	L		M EN			G
VL	R		M EN			G

Additional Accessories:
 Shelltread

Caucasian	C
Tan	T
Brown	B
Jet Black	J

25-30cm

Narrow	N
Wide 25-30cm only	W

College Park suggests the Wide option for male patients.

1-6

WEIGHT KG	0-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125
SIZE CM	25-26								
LOW-MODERATE IMPACT	1	1	1	2	2	3	4	5	N/A
HIGH IMPACT	1	2	2	3	3	4	5	6	N/A
SIZE CM	27-30								
LOW-MODERATE IMPACT	1	1	1	1	2	2	3	4	5
HIGH IMPACT	1	2	2	2	3	3	4	5	6

Erimed International
 Stensättravägen 9A
 127 39 Skärholmen
 Tel: 08- 449 56 50,
order@erimed.se www.erimed.se