

## Beställningsedel College Park *Velocity*


CUSTOMER #	P.O.#	DATE	
<b>PROSTHETIST INFORMATION</b>			
<b>BILLING</b>		<b>SHIPPING (LEAVE BLANK IF SAME AS BILLING)</b>	
FACILITY/ATTN:		FACILITY/ATTN:	
ADDRESS		ADDRESS	
CITY	STATE/PROV	CITY	STATE/PROV
COUNTRY	POST	COUNTRY	POST
PHONE	FAX	PHONE	FAX
CARRIER* <input type="checkbox"/> UPS <input type="checkbox"/> OTHER		DATE REQUIRED	TIME
		EMAIL	
<b>PATIENT INFORMATION</b>			
PROSTHETIST NAME		REQUISITIONER	
PATIENT ID			
NOTES			

**Art.nr. Right:** 761910-1XX

Transtibial  Transfemoral

**Art.nr. Left:** 761910-0XX

Transtibial  Transfemoral



**velocity** **M3**

PART ID	SIDE	SHELL COLOR	MOUNTING	SIZE	WIDTH	FIRMNESS CATEGORY
VL	L		M	EN		G
VL	R		M	EN		G

Additional Accessories:  
 Shelltread

Caucasian	C	25-30cm	Narrow	N	1-6
Tan	T		Wide	W	
Brown	B		25-30cm only		
Jet Black	J		College Park suggests the Wide option for male patients.		

WEIGHT KG	0-45	46-54	55-63	64-72	73-81	82-91	92-100	101-113	114-125
SIZE CM	25-26								
LOW-MODERATE IMPACT	1	1	1	2	2	3	4	5	N/A
HIGH IMPACT	1	2	2	3	3	4	5	6	N/A
SIZE CM	27-30								
LOW-MODERATE IMPACT	1	1	1	1	2	2	3	4	5
HIGH IMPACT	1	2	2	2	3	3	4	5	6