










Dictus Hip

Product picture	Art.no.	Product	Measure and number of parts
	530500-011	Dictus Hip-Knee DS-50M Complete Medium	Hip belt Medium: 55-80 cm Thigh piece Medium: 30-40 cm* 2-piece rubber band: 50cm
	530500-021	Dictus HIP-Knee DS-50L Complete Large	Hip belt Large: 80-130 cm Thigh piece Large: >40 cm* 2-piece rubber band: 50 cm
	530500-311	Dictus HIP-Foot DS-53M Complete Medium	Hip belt Medium: 55-80 cm Foot part barefoot 2-piece rubber band: 100 cm 2 pair S-hooks and 1-piece Monohook
	530500-321	Dictus HIP-Foot DS-53L Complete Large	Hip belt Large: 80-130 cm Foot part barefoot 2-piece rubber band: 100 cm 2 pair S-hooks and 1-piece Monohook
	530510-011	DS-51M Hip belt Medium	Hip belt Medium: 55-80 cm
	530510-021	DS-51L Hip belt Large	Hip belt Large: 80-130 cm
	530520-011	DS-52M Thigh piece Medium	Thigh piece Medium: 30-40 cm*
	530520-021	DS-52L Thigh piece Large	Thigh piece Large: >40 cm*
	530300-011	DS-03 Barefoot	

*Thigh measurement 10 cm above the patella.

	530530-000	Dictus Hip Strap	
	530100-050	Rubber band for Dictus HIP-Knee	Length: 50 cm/single
	530100-101	Rubber band with knee strap for Dictus HIP-Foot	Length: 100 cm/single
	530100-100	Rubber band	Length: 100 cm/single
	530100-003	S-hooks / pair (steel)	
	530100-001	Plasthooks / pair	
	530100-002	Monohook for Dictus / single	

*Thigh measurement 10 cm above the patella.