



Dictus Flex


Improves the mobility of users with foot drop.

The Dictus Band has helped thousands of users who suffer with foot drop to live a normal everyday life again. The band lifts the foot immediately after toe off, which reduces the risk of tripping on doorsteps, the edges of rugs or cobbles. This enables the user to move around unhindered again, enjoy walks and not have to worry about accidents.

The Dictus Band is a discreet aid that is suitable for use with shoes with laces. It consists of an ankle strap, a rubber band and an attachment that is fitted to the shoe.

The ankle strap is secured easily around the leg with velcro and sits comfortably in place all day. The strap is made from real leather, which becomes softer and more comfortable each time you wear it. The ankle strap is connected to the shoe attachment by a silicone rubber band. Dictus Flex features a silicone rubber band with two integrated plates. The Flex silicone rubber band firmly supports the foot and can be easily tightened and adjusted as required.

Read more and watch an instructional video on our website.



Benefits:

- Gives the user confidence and freedom by reducing the risk of tripping and falling.
- Made from durable leather, which gets softer each time you wear it.
- Flex rubber band that can be tightened and adjusted.
- A discreet aid that is suitable for use with most types of shoes.
- Available for both adults and children.

Dictus Flex

	Article no.	Product		Colour
	530400-031	Dictus DS-04XL	Xlarge >26 cm	Black
	530400-021	Dictus DS-04L	Large 20-27 cm	Black
	530400-011	Dictus DS-04M	Medium 15-20 cm	Black
	530400-000	Dictus DS-04J	Child 12-17 cm	Black
	530300-011	DS-03R	Barefoot right	Black
	530300-111	DS-03L	Barefoot left	Black
	530100-400	Flex-elastic strap J/M 35 cm		Black
	530100-440	Flex-elastic strap L/XL 40 cm		Black
	532100-004	S-hooks, 4-pack		
	530100-001	Plastic hooks/pair		Black
	530100-002	Monohook for Dictus/single		Black