



Dictus Hip

Helps Hip and
Ankle Flexors
with Power.

Dictus Hip is suitable for people with weak hip flexors, such as individuals with Multiple Sclerosis (MS) or that had a Stroke, and who have general weakness in their leg(s).

Dictus Hip consists of a comfortable hip belt and an elastic silicone rubber band that is attached either to the thigh piece (Dictus Hip-Knee) or to the foot piece (Dictus Hip-Foot), depending on the user's needs.

The construction assists with lifting the leg when walking, helps the user's balance and ability to move safely.


The silicone rubber bands assist with flexion of the hip. The level and pull direction of the assistance can be easily adjusted to suit the user's needs. It is easy to alter flexion, adduction and the pull direction of inward rotation. If Dictus Hip-Foot is used, foot inversion can also be corrected.

Benefits:

- Helps the user getting the leg into the swing phase of the step, giving them a better gait.
- Provides direct feedback of results.
- Can be worn over ordinary clothing or under loose clothing.
- Universal and can be used on one or both legs.



Dictus Hip

	Article no.	Product
	530500-011	Dictus Hip-Knee DS-50M Complete Medium (hip belt M, thigh piece M, 2 silicone rubber bands 50 cm)
	530500-021	Dictus Hip-Knee DS-50L Complete Large (hip belt L, thigh piece L, 2 silicone rubber bands 50 cm)
	530500-311	Dictus Hip-Foot DS-53M Complete Medium (hip belt M, foot piece, 2 silicone rubber bands 100 cm, 2 pairs S-hooks, Monohook)
	530500-321	Dictus Hip-Foot DS-53L Complete Large (hip belt L, foot piece, 2 silicone rubber bands 100 cm, 2 pairs S-hooks, Monohook)
	530510-011	Dictus Hip DS-51M Hip belt Medium 55-80 cm Waist measurement
	530510-021	Dictus Hip DS-51L Hip belt Large 80-130 cm Waist measurement
	530520-011	Dictus Hip DS-52M Thigh piece Medium 30-40 cm Thigh measurement 10 cm above the patella
	530520-021	Dictus Hip DS-52L Thigh piece Large 40-50 cm Thigh measurement 10 cm above the patella
	530520-031	Dictus Hip DS-52L Thigh piece XLarge > 50 cm Thigh measurement 10 cm above the patella
	530300-011	Dictus in DS-03, foot piece
	530530-000	Dictus Hip Suspenders
	532100-040	Silicone rubber band 40 cm / single - Dictus Hip-Knee
	532100-050	Silicone rubber band 50 cm / single - Dictus Hip-Knee
	532100-080	Silicone rubber band 80 cm / single - Dictus Hip-Foot
	532100-100	Silicone rubber band 100 cm / single - Dictus Hip-Foot
	532100-120	Silicone rubber band 120 cm / single - Dictus Hip-Foot
	532100-121	Knee strap for Dictus Hip-Foot
	530530-001	Knee cap for Dictus Hip-Foot
	532100-004	S-hooks/4-pack
	530100-001	Plastic hooks/pair
	530100-002	Monohook for Dictus/single