

# Dictus Original

Gives users with foot drop confidence and freedom.

The Dictus Band has been helping users who suffer from foot drop to rediscover an active life for many years now. The band lifts the foot immediately after toe off, meaning doorsteps, the edges of rugs and cobblestones no longer pose an obstacle when walking. Users avoid feeling restricted in their everyday life and no longer have to worry about tripping over.


Dictus Original is a discreet aid that is suitable for both adults and children. The strap is secured easily around the ankle with velcro and sits comfortably in place all day. The ankle strap is made from real leather, which

is soft against the skin and becomes more comfortable each time you wear it. There is a D-ring on the front part of the ankle strap, through which the accompanying silicone rubber band is threaded. This is then secured to the hooks attached through the shoes' lacing eyelets. As the user walks wearing the Dictus Band, the foot is automatically lifted by the rubber band. The Dictus Band is available in black and natural leather and is suitable for use with most types of shoes.




Read more and watch an instructional video on our website.

## Benefits:

- Gives the user confidence and freedom by reducing the risk of tripping and falling.
- Made from durable leather, which gets softer each time you wear it.
- Soft rubber band which flexes the foot.
- A discreet aid that is suitable for use with most types of shoes.
- Available for both adults and children.



# Dictus Original

	Article no.	Product		Colour
	530200-031	<b>Dictus DS-02XL</b>	Xlarge >26 cm	Black
	530200-021	<b>Dictus DS-02L</b>	Large 20-27 cm	Black
	530200-011	<b>Dictus DS-02M</b>	Medium 15-20 cm	Black
	530200-000	<b>Dictus DS-02J</b>	Child 12-17 cm	Black
	530300-011	<b>DS-03R</b>	Barefoot right	Black
	530300-111	<b>DS-03L</b>	Barefoot left	Black
	532100-042	<b>Silicone Rubber band, std length 18 cm, 2-pack</b>		Black
	532100-043	<b>Silicone Rubber band, short 16,5 cm, 2-pack</b>		Black
	532100-004	<b>S-hooks, 4-pack</b>		
	530100-001	<b>Plastic hooks/pair</b>		Black
	530100-002	<b>Monohook for Dictus/single</b>		Black