

<b>Orthotech protection shoe</b>			
<b>Unit</b>	<b>B</b>	<b>C</b>	<b>D</b>
Size	Width	Ball	Heel circ.
36	89	239	328
37	91	244	334
38	93	248	342
39	95	251	350
40	96	256	358
41	98	260	368
42	100	265	374
43	101	268	379
44	102	273	388
45	104	277	395
46	105	281	404
47	107	285	409
48	109	289	414
<b>Orthotech Alu.cap. Measurement taken 25% in front of the metatarsal</b>			
<b>Unit</b>	<b>B</b>	<b>E</b>	
Size	Width	Toe height	
36-38	77	46	
39-41	80	48	
42-43	82	50	
44-45	85	52	
46-48	88	53	