



2022
2023

SPORT SHOES

for orthopedic purposes and individual adaptation

Erimed's vision is to give all people with function variations and disabilities an optimal quality of life.

We develop and market products that contributes to ease and optimize the active everyday life of all disabled individuals with the motto "nothing is impossible".

We offer specialist knowledge, high service and engagement, and we have a broad range of high quality elaborated products to orthopaedic related functions.

In this brochure we present our Sport Shoes which is a part of our Comfort Shoe programme. You find our full range of products on our webpage and in other brochures.

As our product range constantly changes, we can not guarantee that all products and sizes in this brochure corresponds to the current situation.

To have a full update on our whole product range, please visit www.erimed.se or contact us at order@erimed.se alternatively tel +4684495650.

Information about the brands presented in this brochure.

Xelero

Xelero shoes are known for their very stable shoes that enables individual adaptation and the outsole provides very good toe-off effect with rocker outsole.

- A. Foam footbed with open cell structure for transportation of perspiration.
- B. Control plate for stability and reduced pressure.
- C. Propulsion Element provides stability, control and promotes forward motion.
- D. Midsole with gel pads for added shock absorption.
- E. Dynamic metatarsal roll bar designed to reduce fore-foot movement and pressure, and provides torsal support under the foot.

Semirigid rubber outsole and removeable insoles.



since 1995
Cambrian
the foot health company

The stable shoes from Cambrian enables individual adaptation and have;

- Extra volume and forefoot width.
- Good toe-off.
- Anatomic last and semirigid outsole.
- Durable leather and/or mesh.
- Removable insoles.

WOMEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Steadfast Black/Periwinkle

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** Wide
Article no: 806006



Steadfast White Coral

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** Wide
Article no: 806042



Steadfast Fit-Go

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker, semirigid outsole, and Fit-Go closure system.

Size: 37-43 • **Last:** Wide
Article no: 806016



All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

WOMEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Ultra Black Lace

Very comfortable sport shoe in leather, with large volume in forefoot, semirigid outsole and good toe-off.

Size: 37-42,5 • **Last:** B (extra depth)

Article no: 831410



Ultra Genesis Silver/Purple

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 36-42,5 • **Last:** B Cambrian

Article no: 830240



Breeze Black

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 36-42,5 • **Last:** B Cambrian

Article no: 832500

**In stock from
Autumn 2022**



Breeze Grey

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 37-42,5 • **Last:** B Cambrian

Article no: 832510

**In stock from
Autumn 2022**

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

WOMEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Matrix Mesh

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** B

Article no: 805821



Arandi Lace

Stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** D

Article no: 803110



Arandi Velcro

Stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** D

Article no: 803314

Arandi Narrow Velcro

Stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** B

Article no: 804887



Genesis XPS Salmon

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 37-43 • **Last:** B spec

Article no: 802422

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

WOMEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.

**In stock
from Winter
2022/2023**



Steadfast Walker Black

Modern, stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.
Size: 37-44 • Last: B
Article no: 807400

Steadfast Walker Black

Modern, stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.
Size: 37-44 • Last: D
Article no: 809740



Steadfast Walker White

Modern, stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.
Size: 37-44 • Last: D
Article no: 807401

**In stock
from Winter
2022/2023**



Steadfast High Hiker

Stable and comfortable waterproof hiker boot in nubuck with very good toe-off / rocker and rigid outsole.
Size: 37-43 • Last: Medium
Article no: 801252

**In stock from
late Autumn
2022**



Steadfast Trail

Stable and comfortable water-proof trail shoe in mesh fabric with very good toe-off / rocker and rigid outsole.
Size: 37-43 • Last: Medium
Article no: 800446

**In stock from
late Autumn
2022**

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

MEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Steadfast Black

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 41-48 • **Last:** Wide

Article no: 802800



Steadfast Fit-Go

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker, semirigid outsole, and Fit-Go closure system.

Size: 41-48 • **Last:** Wide

Article no: 802860



In stock
from Winter
2022/2023



Steadfast Walker Black

Modern, stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 41-49 • **Last:** D

Article no: 808300

Steadfast Walker Black

Modern, stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 41-49 • **Last:** 4E

Article no: 805830

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

MEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Ultra Genesis Silver/Blue

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 41-47 • **Last:** 2E

Article no: 830142



Bora Grey

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 41-47 • **Last:** D

Article no: 832410

**In stock from
Autumn 2022**



Bora Black

Very comfortable sport shoe in mesh fabric with good volume in forefoot, semirigid outsole and good toe-off.

Size: 41-47 • **Last:** D

Article no: 832400

**In stock from
Autumn 2022**



Ultra Black Lace

Very comfortable sport shoe in leather, with large volume in forefoot, semirigid outsole and good toe-off.

Size: 41-48 • **Last:** 2E (extra depth)

Article no: 831312

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

MEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.



Arandi Lace

Stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 41-48 • **Last:** 4E

Article no: 805110



Arandi Velcro

Stable and comfortable sport shoe in leather with very good toe-off / rocker and semirigid outsole.

Size: 41-48 • **Last:** 4E

Article no: 805114



Genesis XPS Grey

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 41-48 • **Last:** D spec

Article no: 804630



Matrix Mesh

Stable and comfortable sport shoe in mesh fabric with very good toe-off / rocker and semirigid outsole.

Size: 41-48 • **Last:** D

Article no: 805815

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

MEN'S SHOES

All shoes have removable insoles. All products in this brochure are stocked. Please also note that products temporarily can go out of stock. We reserve the right to assortment changes and possible printing errors.

In stock from
late Autumn
2022



Steadfast Trail

Modern, stable and comfortable waterproof trail shoe in mesh fabric with very good toe-off / rocker and rigid outsole.

Size: 41-48 • **Last:** Medium
Article no: 807257

In stock from
late Autumn
2022



Steadfast High Hiker

Stable and comfortable waterproof hiker boot in nubuck with very good toe-off / rocker and rigid outsole.

Size: 41-48 • **Last:** Medium
Article no: 803000

Steadfast High Hiker

Stable and comfortable waterproof hiker boot in nubuck with very good toe-off / rocker and rigid outsole.

Size: 41-48 • **Last:** Extra Wide
Article no: 807300

All models are different – thus – Please always look at the last measurement tables for reference, pages 10-11.

Arandi Women – D last

Unit Size	A Length	B Width	C Ball
37	242	86	230
38	250	88	237
39	259	90	243
40	263	91	246
41	271	93	252
42	276	94	256
43	284	96	262

Arandi Men – 4E last

Unit Size	A Length	B Width	C Ball
41	270	99	272
42	274	100	275
43	283	101	281
44	287	103	285
45	296	105	291
46	300	106	294
47	308	108	300
48	317	110	307

Matrix, Arandi Narrow Women – B last

Unit Size	A Length	B Width	C Ball
37	242	83	221
38	250	85	227
39	259	87	233
40	263	88	237
41	271	90	243
42	276	91	246
43	284	93	252

Matrix Men – D last

Unit Size	A Length	B Width	C Ball
41	270	93	254
42	274	94	257
43	283	96	263
44	287	97	266
45	296	99	273
46	300	100	276
47	308	103	282
48	317	105	289

Genesis XPS Women – B-spec last

Unit Size	A Length	B Width	C Ball
37	241	84	230
38	249	87	237
39	258	89	243
40	262	90	246
41	270	92	253
42	275	93	256
43	283	95	262

Genesis XPS Men – D-spec last

Unit Size	A Length	B Width	C Ball
41	275	93	256
42	279	94	259
43	287	96	265
44	292	97	269
45	300	99	275
46	305	100	278
47	313	102	284
48	321	104	291

Steadfast Women – Medium & B last

Unit Size	A Length	B Width	C Ball
37	234	84	222
38	243	86	229
39	251	88	235
40	255	89	238
41	264	91	245
42	268	92	248
43	276	94	254
44	285	96	260

Steadfast Women – Wide and D last

Unit Size	A Length	B Width	C Ball
37	234	85	231
38	243	87	238
39	251	89	244
40	255	90	247
41	264	92	254
42	268	93	257
43	276	95	263
44	285	97	269

Steadfast Men – Medium & D last

Unit Size	A Length	B Width	C Ball
41	274	94	254
42	278	95	257
43	286	97	263
44	290	98	266
45	299	100	273
46	303	101	276
47	312	103	282
48	320	106	289
49	329	108	295

Steadfast Men – Wide & 2E last

Unit Size	A Length	B Width	C Ball
41	275	95	271
42	279	96	274
43	287	98	280
44	291	99	283
45	300	101	290
46	304	102	293
47	313	104	299
48	321	107	306
49	330	109	312

Steadfast Men – Extra Wide & 4E last

Unit Size	A Length	B Width	C Ball
41	275	96	275
42	279	97	278
43	288	99	284
44	292	100	287
45	300	102	294
46	305	103	297
47	313	105	303
48	322	108	310
49	330	110	316

Ultra Black – 2E last Men

Unit Size	A Length	B Width	C Ball
41	275	100	264
42	279	101	268
43	287	103	274
44	292	104	277
45	296	105	280
46	304	107	287
47	313	109	293
48	321	111	299

Ultra Black – B Women

Unit Size	A Length	B Width	C Ball
37	250	89	234
38	254	90	237
39	258	91	240
40	267	93	247
41	275	95	253
42.5	284	97	259

Ultra Genesis Men – 2E last


Unit Size	A Length	B Width	C Ball
41	275	100	258
42	279	101	262
43	287	103	268
44	292	104	271
45	296	105	274
46	304	107	281
47	313	109	287

Breeze, Ultra Genesis Women – B last

Unit Size	A Length	B Width	C Ball
36	242	87	222
37	250	89	228
38	254	90	232
39	258	91	235
40	267	93	241
41	275	95	248
42.5	284	97	254

Bora Men – D last

Unit Size	A Length	B Width	C Ball
41	275	98	249
42	279	99	253
43	287	101	259
44	292	102	262
45	296	103	265
46	304	105	272
47	313	107	278



Measures taken on last

**We want to give all people
with function variations
and disabilities an optimal
quality of life.**

Read more about us and our
products at www.erimed.se



www.erimed.se



Ellipsvägen 11

141 75 Kungens Kurva

Phone: +46 (0)8 - 449 56 50

Email: info@erimed.se

Graphic Design: ClearDesign.se