

Orthotech

Orthotech Stabil Safety				
Unit	A	B	C	D
Size	Length	Width	Ball	Heel
36	234	89	239	328
37	243	91	244	334
38	251	93	248	342
39	259	95	251	350
40	264	96	256	358
41	268	98	260	368
42	276	100	265	374
43	285	101	268	379
44	290	102	273	388
45	296	104	277	395
46	302	105	281	404
47	310	107	285	409
48	318	109	289	414

Orthotech Stabil, Reha, Vario Incl. thick resp. thin insole							
Unit	Referens	A	B	C- thick	C- thin	D- thick	D- thin
UK stl	EU stl	Length	Width	Ball	Ball	Heel	Heel
3	35,5	234	83	227	230	326	329
4	37	243	85	233	236	337	340
5	38	251	87	238	241	343	346
6	39	259	89	243	246	355	358
7	40,5	268	91	248	251	368	371
8	42	276	93	254	257	374	377
9	43	285	95	259	262	381	384
10	44,5	293	97	264	267	389	392
11	46	302	99	269	272	399	402
12	47	310	101	273	276	411	414
13	48	318	103	279	282	416	419
14	50	327	105	288	291	428	431